

Listen to what some of our Satisfied Customers have to say about the therapeutic value of the VacuPractor!

Hi Paul,

My back went "out" a few weeks ago and I have to say that the vacupractor has helped more than I can believe.



The neat thing is as you get more comfortable using it you can start doing little things that make it even more effective. I have gotten to where after a few minutes I flex my lower back just a tiny bit, it increases the suction and pulls harder; after about 45 seconds I relax and the relief is immediate.

Ken C.

Hey Paul!

Wow! Thank you for The Vacu Thing! Amazing results in a very short time. After about 5 uses, I was noticing a remarkable decrease in lower back pain. Also, while riding my bicycles I have experienced a significant increase in power! Most noticeable while climbing the steep hill's here in the Seattle area, which range up to 20%. Thank You very much. I am standing straighter and taller. I'm pain free. And I can ride faster and further. Remarkable. Wow.

Scott F.

To learn more about how the VacuPractor can improve the quality of your life,

Call: 800-676-9320
or Visit:


VacuPractor.com



- A patented medical device registered with the FDA.
- Complimentary to Chiropractic care.
- Increases the clearances between the vertebrae.
- Can eliminate the need for muscle-relaxers and surgery.
- Can reverse or slow the onset of lower spinal problems.



15127 NE 24th Street #185
Redmond, WA 98052



Eliminate back pain, and feel like a new person, in just 5 to 15 minutes a day...

GUARANTEED!

WITH





The VacuPractor™ – How It Works!



1. Air Pocket Sealed

The VacuPractor™ is used at home on a flat firm surface.



2. Trapped Air Expelled

1. Lie down on the VacuPractor™ after wetting it with a damp cloth or spray bottle. A pocket of air will be trapped behind the curvature of the lower spine and the VacuPractor.



3. Vacuum Pulls Spine Flat

2. Raise your knees to chest, one time to expel the trapped air and gently align the spine.



4. Pain Relieved!

3. Lower your feet to the floor, a vacuum is formed that pulls straight down on the curvature of the lower back.

4. Straighten your legs to lengthen the spine and relieve pressure on the discs.

Relax on the VacuPractor™ for 5 to 15 minutes, enjoying the gentle vacuum pressure. You'll feel the back pain slip away.



So Simple, Yet So Effective!

“From the first time I used the VacuPractor, I knew the treatment of back pain had been changed forever. Now relief is at your finger tips with this revolutionary new technology that my patients just can't seem to get enough of.”


- Dr. Eric Hansen (Doctor of Chiropractic)

Benefits of The VacuPractor™

1. Relieves pain & stops the inflammation cascade.
2. Increases the clearances between the vertebrae for better nerve function.
3. Stretches the muscles, reducing or eliminating the need for muscle relaxing medications or expensive and painful surgery!
4. Reverses or slows the onset of lower spinal problems.
5. **Get back to living a normal, pain-free life!**



Michael George, Trainer to the Stars, Personal Health Coach & Fitness Expert



Michael George, Trainer to the Stars, Personal Health Coach, Fitness Expert and Author of “Body Express Makeover.”

Michael is well known for his fitness training work with such celebrities as actors Dennis Quaid, Reese Witherspoon, Meg Ryan, Christian Slater, Tobey Maquire, record producer and rapper P. Diddy, musician Slash, and many others. In addition, Michael is the author of “Body Express Makeover” and has been featured on CNN, Fox, Healthy Lifestyle and The TV Guide.

After personally trying The VacuPractor Michael had this to say about it:

“I just wanted to let you know I received the VacuPractor. I tried it for a couple of days and I was really amazed at how such a simple concept could help decompress my back so well. I could really feel the vacuum compression and my back feels as good as ever. This is a great product that will help a lot of people. I am thoroughly impressed and I'm all in!”